



Last modified: March 2019  
Last reviewed: March 2019

## Protocol: Isolation of genomic DNA with NucleoSpin Blood L – Midi Kit

**Purpose:** This protocol is used to extract gDNA from frozen cell pellets using the NucleoSpin Blood L Midi Kit. The manufacturer's protocol recommends this kit for cell pellets between  $5 \times 10^6$  –  $2 \times 10^7$  cells. It is strongly encouraged to first test this kit on not-infected cells to ensure this same range applies to the cells of interest. gDNA recovery from the column can be estimated using an average value of ~6.6 pg gDNA per cell.

### Materials:

- NucleoSpin Blood L – Midi Kit (Clontech Cat. #740954.20)
- RNase A (Clontech, Cat. #740505), lyophilized, dissolve to 20 mg/ml in H<sub>2</sub>O; store at -20°C
- Absolute ethanol
- Centrifuge with a swing-out buckets rotor capable of reaching 3500 x g is required

### Getting started:

- Set a water bath to 56°C.
- Preheat Elution Buffer BE on top of a heat block set to 70°C.
- Equilibrate cell pellets to room temperature (15–25°C).
- Prepare Buffer BQ2 by adding ethanol and prepare Proteinase K by resuspending in Proteinase Buffer according to the kit insert.
- Label one 15 mL centrifuge tube for every column that will be used. *These tubes do not have to be the specific collection tubes provided in the NucleoSpin kit. Any 15 mL conical (e.g. Falcon™ 15mL Conical Centrifuge Tubes) can be used.*

### Step 1. Lyse cell pellet

1. Add 1 mL PBS and 150  $\mu$ L of Proteinase K to each conical.
  - 1 mL PBS assumes the cell pellet volume is 1 mL, bringing the final volume to approximately 2 mL.
2. Add cell pellet and mix well by vortexing.
3. Add 4.1  $\mu$ L RNase A (20 mg/mL) to the 15 mL conical and vortex.
  - If the cell pellet is still clumpy, incubate at 56°C for 5-20 min and vortex the tube several times to further break up the cells.
4. Add 2 mL of Buffer BQ1 and vortex the mixture vigorously for 10s.

- Vigorous mixing is important for DNA yield and purity.
- 5. Incubate samples at 56°C in the water bath for 15-45 minutes, vortexing every 10 minutes until solution is clear.
- 6. Let the lysate cool down to room temperature before proceeding with ethanol precipitation.
  - Make sure that the lysate has cooled down to room temperature (about 5 min) before loading it onto the columns. Loading of hot lysate may result in reduced yields due to evaporation of the ethanol.

## Step 2. Precipitate DNA

- 7. Add 2 mL absolute ethanol to each sample and *immediately* mix by inverting the tube 10 times.
  - High local ethanol concentration must be avoided by immediate mixing.

## Step 3. Bind DNA

- 8. For each prep, take one NucleoSpin Blood L Midi column placed in the provided conical and load 3 mL of lysate. Avoid drops on the rim of the column. Close the tubes and centrifuge 3 min at 3250 x g. Discard flow through.
  - The lysate will usually start to flow through the column even before centrifugation. This will not adversely affect DNA yield or purity.
  - Keep NucleoSpin Blood L Midi columns in an upright position as liquid may pass through the ventilation slots on the rim of the column even if the caps are closed.
- 9. Load the remaining 3 mL of lysate to the respective NucleoSpin Blood L Midi columns. Avoid drops on the rim of the columns. Close the tubes and centrifuge 5 min at 3250 x g.
  - Remove the collection tubes with the columns carefully from the rotor and avoid having the flow-through come in contact with the column outlets.
- 10. Discard the flow-through and place the column back into the collection tube.

## Step 4. Wash silica membrane

- 11. **1st wash:** add 2 mL of Buffer BQ2 to the NucleoSpin Blood L Midi columns. Centrifuge 2 min at 3250 x g. It is not necessary to discard the flow-through after the first washing step.
- 12. **2nd wash:** add 2 mL of Buffer BQ2. Centrifuge 10 min at 3250 x g. Remove the columns carefully from the rotor to avoid having the flow-through contact the column outlets. Discard flow through.
  - A longer centrifugation during this second washing step insures that residual ethanol from Buffer BQ2 is removed from the silica membrane of the column so no further drying step is necessary.

## Step 5. Elute DNA

See the Practical Considerations Appendix for alternative elution conditions that may improve gDNA yields.

13. Insert the column into a new 15 mL collection tube and apply 200 uL of preheated (70°C) Elution Buffer BE directly to the center of the silica membrane.
15. Incubate overnight at 4°C, or at room temperature for 5 minutes.
16. Centrifuge at 3250 x g for 2 min to elute.

## Step 6. Dilute and store gDNA

17. Measure the concentration of gDNA (for example, using NanoDrop or Qubit). If concentration is higher than 200 ng/uL, dilute gDNA with Elution Buffer BE.
  - We recommend calculating the dilutions based on a precisely measured elution final volume because there is always some loss due to the column.
18. Store at 4°C for short term, or -20°C for long term.
19. Next, please see: **“Guidelines for gDNA sample submission for PCR and sequencing.”**

## **Practical Considerations Appendix**

### **RNase:**

The original NucleoSpin Blood L protocol does not call for the addition of RNase. However, the presence of RNA contamination in gDNA samples can skew quantification when using absorbance-based methods for determining sample concentrations. Furthermore, large amounts of RNA contamination can have a negative impact on subsequent PCR steps. For these reasons we highly recommend using RNase when performing large scale gDNA isolations.

### **Elution Volume:**

The Genomic DNA from Blood User Manual lists the following elution procedures for highest yield and/or concentration:

- High yield: Perform two elution steps with the volume indicated in the individual protocol. About 90–100% of bound nucleic acid can be eluted
- High concentration: Perform one elution step with 60% of the volume indicated in the individual protocol. Concentration of DNA will be higher than with standard elution. Maximal yield of bound nucleic acid is about 80%.
- High yield and high concentration: Apply half the volume of elution buffer as indicated in the individual protocol, incubate for 3 min and centrifuge. Apply a second aliquot of elution buffer, incubate and centrifuge again. Thus, about 85–100% of bound nucleic acid is eluted in the standard elution volume at a high concentration.