

The RNAi Consortium

TRC Laboratory Protocols

Protocol Title: **Clone oligos into pLKO vectors for shRNA Constructs**

Current Revision Date: **2/12/2013**

RNAi Platform, Broad Institute, trc_info@broadinstitute.org

Brief Description:

This protocol is adapted from TRC high Throughput shRNA Cloning Protocol to allow users to clone synthetic oligos into pLKO as RNAi constructs (on a small scale).

Table of Contents:

- I. Materials
- II. Instruction
- Revision Notes

These procedures should be carried out in accordance with biosafety requirements of the host institution.

Protocol:

I. Materials

- TB media (terrific broth; American Bioanalytical,)
- Carbenicillin (Sigma)
- T4 DNA ligase, 400units/ul (NEB)
- 10X ligase buffer (NEB)
- Restriction enzymes (NEB)
- Qiaquick gel extraction kits (Qiagen)
- DH5 α Max Efficiency Cells (Invitrogen, or Gene Choice)
- SOC Media (Invitrogen)
- LB agar plates (with appropriate antibiotic)
- Big dyes (ABI)
- Stuffed-pLKO vector, pLKO_TRC001 or pLKO-TRC005 (TRC)
- Hairpin oligonucleotides (IDT)
- 10X Buffer 2 (NEB)

II. Instructions

1. Vector preparation

Circular stuffed pLKO vector is digested by AgeI and EcoRI:

- In a 1.5ml eppendorf tube, add:
 - 20ug of circular vector, stuffed-pLKO
 - 20ul NEB buffer I
 - Appropriate amount of H₂O
 - 6ul AgeI
 - 6ul EcoRI
 -
 - 200ul total reaction
- Set up two controls as well; controls will be loaded in the purification gel:

The RNAi Consortium (TRC)

Broad Institute

A. AgeI control:

1 μ g circular vector
2ul NEB buffer 1
Appropriate amount of H₂O
1ul AgeI

20ul

B. EcoRI control:

1 μ g circular vector
2ul NEB buffer 1
Appropriate amount of H₂O
1ul EcoRI

20ul

Incubate at 37°C for 4-6 hours.

- Fragments are purified by a 1% agarose gel (120v for 17 hours). Add 40ul 6x loading dye into the 200ul AgeI/EcoRI digestion and load 15ul sample in each well. Restriction controls, intact plasmid, and linear size markers are loaded in a single well each.
- The 7kb open-vector band is excised and the vector is extracted using Qiaquick Gel Extraction kit (Qiagen). The extracted DNA is suspended into a final volume of 2000ul of TE buffer (about 6 ng per ul). The opened vector is then stored frozen.

Note: When stuffed-pLKO is used to prep the linear vector, its 1.9kb stuffer facilitates the high-grade purification by separating any incompletely digested vectors from double-digestion products. As a result, the cloning background (ie. empty vector) is no longer a concern.

2. Hairpin oligonucleotides features

Oligos should be designed to have following configuration. The forward and reverse strands pair with each other and result in 2 overhangs that can be ligated to the processed pLKO vector. To avoid insert chimerism, we do not recommend adding phosphates to oligo strands.

Forward oligo: 5' CCGG---21 bp Sense---CTCGAG---21 bp Antisense---TTTTTG3'
Reverse oligo: 5' AATTCAAAA---21 bp Sense---CTCGAG---21 bp Antisense---3'

3. Anneal oligos

- Re-suspend oligos in ddH₂O (100uM)
- In an eppendorf tube, add:

The RNAi Consortium (TRC)

Broad Institute

1.5 µL of S oligo (100uM)
1.5 µL of AS oligo (100uM)
5 µL of 10x NEB buffer 2
42 µL ddH₂O

50ul total

- Incubate for 4 min at 95°C.
- Incubate for 10 min in 1L beaker filled with 70°C water and allow the water in the beaker to cool to RT (*Note: This will take a couple of hours, but it is important that you do this slowly! One may set up a thermal cycler to achieve the same.*)

4. Ligation

- Each ligation contains the following to be added in the order listed:

Note: Users may make a serial dilution e.g. 1ul straight, 1:10, 1:100, 1:1000, to capture an appropriate stoichiometry (between oligo and vector) in the ligation reaction.

1ul of annealed oligo pair
3ul of prepped open vector (~20ng)
2ul 10x ligase buffer (NEB)
13ul of H₂O
1ul T4 ligase (NEB)

20ul total

- Ligation proceeds at 16°C for 3-4 hours.

5. Transformation/plating/picking colonies

- 2ul ligation reaction mix is added to 25ul competent cells (DH5α). The mix is placed on ice for 30min.
- Heat shock at 42°C for 90 seconds, followed by incubation on ice for 2min.
- The transformed cells are recovered with 400ul SOC at 37°C for 1 hour.
- ~150ul transformants are plated onto agar dishes containing 100ug/ml carbenicillin.
- Incubate agar dishes upside down at 37°C overnight.
- Pick 2-4 colonies each into 1-2 ml TB + 100ug/ml carbenicillin.
- Cells grow at 37°C (300rpm) for no more than 17 hours.
- Glycerol stock can be made by mixing 150ul glycerol with 850ul culture. Stored at -80°C.

6. Hairpin sequencing protocol

- A. Primer sequence (for both pLKO_TRC001 and pLKO-TRC005):

5'-GAT ACA AGG CTG TTA GAG AGA TAA TT-3'

- B. DNA prep: DNA is prepped by TempliPhi from glycerol stock, or by alkaline lysis prep from overnight cultures. Adjust the concentration to 40ng/ul.

The RNAi Consortium (TRC)

Broad Institute

C. Sequencing Reaction: Each sequencing reaction contains ~50ng DNA template. We use 4:1 ratio of bigDye3.1:dGTP bigDye to read through the hairpin structures.

D. Cycled as following:

96°C/45" -- [96°C/15" - 50°C/15" - 60°C/3' 30"] x 34 cycles -- hold at 4°C

E. Sequencing reactions are cleaned up by ethanol precipitation. The dried fragments are eluted into 10ul of 0.1 mM EDTA, detected with ABI 3730xl, with POP7, 36cm array (Injection voltage 1.2, Injection time 15 sec, Run time 2450 sec, Oven temperature 60, Run voltage 8.5).

Revision Notes: