

Last modified: August 2018

Last revised: August 2018

## Protocol: Assaying dCas9 Activity for CRISPRa screening with Adherent Cells

After having established a dCas9-expressing cell line, this assay is to determine the level of dCas9 activity. A sgRNA for CD4 and a sgRNA for CD45 is introduced using a lentivirus (XPR\_502). When CD4 and CD45 sgRNA's are introduced into a dCas9-expressing cell line, expression of the corresponding target should increase in comparison to the baseline expression.

Before starting, it is recommended you complete the following steps:

- Puromycin or Blasticidin dosing (see **Protocol: Puromycin and Blasticidin titration**)

Materials and Reagents required:


9 million dCas9 (pXPR\_BRD109) cells  
sgCD4 in XPR\_502 lentivirus (1x1 mL aliquot)\*  
sgCD45 in XPR\_502 lentivirus (1x1 mL aliquot)\*  
Polybrene  
Puromycin  
96-v clear bottom plate (Corning, #3357)  
6-well plates  
12-well plates  
APC anti-human CD4 antibody (Biolegend, #357408)  
FITC anti-human CD45 antibody (Biolegend, #368508)  
Flow buffer (PBS, 2% FBS, 5 $\mu$ M EDTA)  
Flow cytometer (Accuri)

\*The pXPR\_BRD109, sgCD4, and sgCD45 lentiviruses are available from GPP ([gpp-reagents@broadinstitute.org](mailto:gpp-reagents@broadinstitute.org)).

## DAY 1

1. Trypsinize and count cells. Resuspend 9 million cells in a final volume of 9 mL of media.
2. Add 12 uL of 1000x polybrene.
3. Add 1.5 mL of cell suspension per well into 5 wells of a 12-well plate.
4. Add virus to wells.
  - a. Add 500 uL of media to 1 well [no-infection control (NIC)].
  - b. Add 100 uL of sgCD4 virus and 400 uL of media to 1 well.
  - c. Add 500 uL of sgCD4 virus to 1 well.
  - d. Add 100 uL of sgCD45 virus and 400 uL of media to 1 well.
  - e. Add 500 uL of sgCD45 virus to 1 well.

Example of infection plate:




5. Spin plate for 2 hours at 1,000 x g at 30°C. After spin is complete, add dropwise 2 mL of media per well (for a total volume of 4 mL per well) and return plate to incubator.
6. 4-6 hours after the spin is complete, transfer the cells from each well into a T75 or 6-well dish (depending on optimal seeding density).

## DAY 3

1. Trypsinize and count cells.
2. Seed an in-line titer into 6-well plates. For each cell population, seed 2 wells with 1.5E5 cells per well (for a total of 10 wells). Add puromycin to half of the wells (one for each population).

Example of in-line titer plate:



3. For the NIC cells, seed remaining cells in media without puromycin in a T75 or 10 cm dish at standard confluence.
4. For the sgCD4- and sgCD45-infected cells, seed remaining cells in puromycin-containing media in a T75 or 10 cm dish at approximately 50% confluence (to account for selection).

#### DAY 4-8

1. Passage and/or maintain the in-line titer and flasks. Make sure that the - puro wells in the in-line titer do not reach confluence.

#### DAY 8

##### Read out in-line titer

1. Visually confirm that all cells in the + puro NIC well are dead. If they are not, continue passaging titer and puromycin-selection flasks until they are.
  - a. This may also indicate that the puromycin dose is too low and should be re-determined.
2. Trypsinize and count the cells in each well. Obtain the infection efficiency for each virus volume by dividing the number of cells in the + puro well by the number of cells in the - puro well. For each virus, select the volume that gave the closest to 15-40% infection efficiency to stain and flow in the following step.

##### Stain and flow

1. Trypsinize and count cells from the unselected NIC flask and the chosen virus volumes of sgCD4-infected and sgCD45-infected cells.
2. In a 96 well V-bottom plate, seed approximately  $2.5 \times 10^5$  cells/well of the following:
  - a. 2 wells of NIC cells (1 unstained control and 1 stained control)
  - b. 1 well of sgCD4-infected cells

- c. 1 well of sgCD45-infected cells
3. Spin plate (5 min @ 1,000 x g) and remove supernatant with aspirating wand.
4. Turn off hood and room light. Resuspend each well in 90 uL of flow buffer.
  - a. Add 10 uL of flow buffer to the NIC unstained control well.
  - b. Add 5 uL of anti-CD4 antibody and 5 uL of anti-CD45 antibody to each remaining well, for a total volume of 100 uL per well.
  - c. Resuspend each well.
5. Cover the plate with foil and incubate on ice for 20-30 min.
6. Spin (5 min @ 1,000 x g) and remove supernatant with aspirating wand.
7. Add 200 uL flow buffer/well and resuspend.
8. Repeat steps 6 – 7 two more times.
9. Flow plate using a flow cytometer. If using the Accuri use: FL1-A for FITC-anti-CD45 and FL4-A for APC-anti-CD4.
  - a. Gate live cells using FSC/SSC.
  - b. Use the NIC unstained control well to set gates on FL1-A and FL4-A. CD4 and CD45 activation should be assessed with both the % positive cells and the mean fluorescence.