

GENETIC PERTURBATION PLATFORM

Last modified: January 2018 Last reviewed: February 2019

Protocol: Lentivirus concentration protocol

Brief Description:

This protocol describes a method to concentrate lentiviral preps using spin columns. A range of 50 to 200-fold concentration and a volume of 300-500 uL per 60 mL of unconcentrated viral supernatant have been obtained.

Materials & Instrumentation:

- Centrifuge with swinging buckets able to reach 4,000 g that can maintain a temperature of $30^{\circ}C$
- Centrifuge buckets able to hold filters (Ex: Beckman Coulter Inc. #349946)
- Centrifuge bucket snap-on lids (for safety) (Ex: VWR #5000226124)
- Centricon Plus-70 Centrifugal Filter, Ultracel-PL Membrane, 3kDa (Sigma Aldrich #UFC700308)
- Autoclaved Milli-Q water
- 70% ethanol
- 10% bleach
- Lentivirus stock, harvested out of flasks and clarified (1 min spin at 1000 rpm, remove viral supernatant leaving small 293T pellet) right before adding to the columns





Sterilization of filters:

- 1. Snap the filter cups into the collection cups, add 30 mL of 70% ethanol to each filter cup, seal with the filter lids and let sit for 5 min. Centrifuge at 4,000 g for 25 min at 30°C.
- 2. Invert the filter cups in the concentrate cups and centrifuge at 4,000 g for 5 min at 30°C to remove any excess ethanol.
- 3. Reassemble the columns and add 20 mL of autoclaved Milli-Q water to each filter cup. Seal with the filter lids and centrifuge at 4,000 g for 10 min at 30°C to wash the filters.
- 4. Invert the filter cups in the concentrate cups and centrifuge at 4,000 g for 5 min at 30°C to remove any excess water.
- 5. Repeat steps 3-4.

Virus concentration:

- 6. Label the side of the filter cups with the name/condition of the virus.
- 7. Reassemble the columns and add 60 mL of autoclaved Milli-Q water to each filter cup. Seal with the filter lids and centrifuge at 4,000 g for 1 hr at 30°C to hydrate the filters.
 - During this spin harvest the virus from the flasks:
 - Pipette off virus from the flask(s) into conical tube(s)
 - Clarify virus by centrifuging at <u>1000 rpm for 1 min</u>
 - Pipette off viral supernatant, being careful to not disturb any cell pellet
 - Filter the virus if necessary with 0.45 uM filter.
- 8. Invert the filter cups in the concentrate cups and centrifuge at 4,000 g for 5 min at 30°C to remove any excess water.
- 9. Sterilize the lids of the filters in a bucket of 70% ethanol and let dry.
- 10. Reassemble the columns and add up to 60 mL of viral supernatant to each filter cup. Seal with the filter lids and centrifuge at 4,000 g for 1 hr at 30°C.
- 11. Sterilize concentrate cups in a bucket of 70% ethanol and let dry.

- 12. Invert the filter cups into the concentrate cups and centrifuge at <u>1,000 g for 2 min</u> at 30°C to recover the concentrated virus.
- 13. If concentrating multiple columns of the same virus, combine all of the concentrated virus together, then aliquot.
- 14. Bleach all filters cups, collection cups, concentrate cups, filter lids, conicals, filters and flasks for at least 20 min.