

Last modified: June 2016  
 Last reviewed: June 2016

## Protocol: No-Spin Infection Protocol for adherent cell lines

While the Spinfection protocol (ie includes a centrifugation step after the addition of virus) uses lower virus volumes per number of cells infected, this alternative protocol reduces handling time, necessity of centrifugation and over-confluence issues present in the Spinfection protocol.

The GFP-expressing, puromycin/blastocidin-resistant pRosetta virus is available from GPP or the plasmid can be obtained from Addgene (plasmid # 59700).

### Procedure

#### 1. For a New Cell Line:

##### Step 1. Determine Number of Cells to Seed

# of cells/well *	Total Vol in each well (mL)
2.00E+04	1.1
5.00E+04	
1.00E+05	
3.00E+05	
6.00E+05	
1.50E+06	

\*This range of values is for cell lines without especially large or small size and doubling times of ~1 day

**Day 0, pm: Seed cells in 6 wells of a 6-well plate**

**Day 1, am (17hr post-seed)**

Aspirate viral media and add 2.75 mL media per well

**Day 2, pm (48hr post-seed)**

Check confluence of all wells.

Use condition that is 90-95% confluent.

# of cells seeded per well on Day 0	% confluence 48hr post-seed
1.00E+05	
2.50E+05	
5.00E+05	
7.50E+05	
1.00E+06	
2.00E+06	

Examples of appropriate cell numbers at time of seeding:


- A375: ~ 2.5e5 cells
- HT29: ~ 1.0e6 cells
- 293T: ~ 2.0e5 cells

## Step 2. Determine optimal Polybrene (PB) concentration

### Day 0, pm: Seed cells

- 1) Create 1mL of 0.05 ug/uL PB stock
- 2) Create a cell + media master mix:
  - Master mix is made in 20% excess:  $1.2 \times 12 \text{ ml} = 14.4 \text{ ml}$
  - Total volume:  $14.4 \times 0.800 \mu\text{L} = 11.5 \text{ ml}$
  - Cells:  $14.4 \times (\# \text{ of cells from Part 1 Step 1})$
- 3) Seed 12 wells in two 6-well plates with 800uL of cell + media master mix.

PB Condition	Media + Cells (uL)	0.05 ug/uL PB Stock (uL)	Additional Media (uL)	Virus (uL)
0	800	0	140	160
			300	-----
0.5	800	11	129	160
			289	-----
1	800	22	118	160
			278	-----
2	800	44	96	160
			256	-----
4	800	88	52	160
			212	-----
6	800	132	8	160
			168	-----



- 4) Add appropriate volume of PB stock, additional media and virus

**Day 1, am (17hr post-seed)**

Aspirate media from all wells and add 2.75mL new media to each well

**Day 2, pm (48hr post-seed)**

- 1) Determine PB toxicity & infection efficiency
  - Count cells to determine PB toxicity
  - Flow for % GFP+ (or other easily measured marker) to determine efficiency
- 2) Looking at the infection efficiencies and cell counts, determine the PB concentration for your cell line that best balances increased infection efficiency with a toxicity level that you are willing to accept.
- 3) If no PB concentration works, repeat using higher virus volume. Reduce the amount of media used so that the total volume is 1.1mL per well.

**2) For a new virus - Determine virus volume to use**

**Step 1. Infect cells with different virus vol.**

**Day 0, pm: Seed cells**

- 1) Create 2mL of PB Stock at 10X of PB concentration selected in **Part 1 Step 2.**
- 2) Create a cell + media master mix using:
  - Master mix is made in 20% excess:  $1.2 \times 12 \text{ ml} = 14.4 \text{ ml}$
  - Total volume:  $14.4 \times 0.6 \text{ ml} = 8.65 \text{ ml}$
  - Cells:  $14.4 \times (\# \text{ of cells from Part 1 Step 1})$
  - PB: 1.6 ml of 10X stock
- 3) Seed 12 wells of 2 x 6 well-plates with 600uL cells-PB-media master mix
- 4) Add additional media so that the final volume after virus addition will be 1.1mL
- 5) Add Virus: 0 (NIC), 100, 200, 300, 400 and 500uL

	Cells-PB-Media Master Mix per well (uL)	Additional Media (uL)	Virus (uL)
No Infection Control (NIC)	600	500	0
100uL Virus		400	100
200uL Virus		300	200
300uL Virus		200	300
400uL Virus		100	400
500uL Virus		0	500

### Day 1, am (17hr post-seed)

Aspirate media from all wells and add 2.75mL new media.


### Day 2, pm (48hr post-seed)

- **If virus expresses GFP or other easily visualizable marker:** flow for %GFP+ or marker
- **If virus confers drug (e.g. puromycin) resistance:** make a titer plate and drug-select infected cells (Explained in **Part 2 Step 2**)

## Step 2. Making a titer plate to assess infection efficiency by drug selection (e.g. Puromycin)

### Day 2, pm (48hr post-seed/end of infection)

- 1) Passage all 6-well plate conditions.
- 2) Seed  $\frac{1}{2}$  number of cells determined in **Part 1 Step 1** in 5mL total vol. in 2 x 6-well plates using 2 wells per virus volume
  - For the 1<sup>st</sup> well of each 2-well set, select with drug (e.g. Puromycin).
  - For the 2<sup>nd</sup> well of each set, don't add drug selection.



### Day 4-6

- 1) Check confluence of wells in the 6-well plates:  
If No Infection Control well *with* Drug Selection is not dead and some wells are approaching confluence, passage cells.
- 2) Titer is ready to read when No Infection Control *with* Drug Selection is dead (*ensure selection is complete*).
- 3) When ready to read, count cells in each well.
- 4) For each 2-well set, determine infection efficiency:  
$$\frac{(\text{\#cells in well *with* selection})}{(\text{\#cells in well *without* selection})} \times 100\% = \text{Infection Efficiency}$$

- 5) Pick a virus volume that is the closest to the efficiency desired.

### 3) Repeating and scaling up using titered virus and previously used cell line

#### Step 1. Infection and Titer Guide

##### **Day 0, pm: Seed Cells**

In T175s, seed  $18.191 \times (\# \text{ cells from Part 1 Step 1})$

- Scale up screen arms

		Tissue Culture Format	6-well	T175
		Surface Area mm <sup>2</sup>	962	17500
Day 0, pm:	Seeding & Infection Start	Time (hr)	0	
		Seed	Part 1 Step 1 value	Part 1 Step 1 value * 18.191
		PB conc. (ug/mL)	Part 1 Step 2 value	
		Virus vol. (mL)	Part 2 Step 1 value	Part 2 Step 1 value * 18.191
		Total infection/Seeding vol. (mL)	1.1	20

18.191 scaling factor based on the surface area of a T175 vs. 6 well:  $\left(\frac{17500\text{mm}^2}{962\text{mm}^2} = 18.191\right)$

##### **Day 1, am (17hr post-seed)**

Aspirate viral media and add 2.75mL media to the 6-well and 50mL to the T175s

		Tissue Culture Format	6 Well	T175
		Surface Area mm <sup>2</sup>	962	17500
Day 1, am	Change Media	Time (hr)	17	
		Media (mL)	2.75	50

##### **Day 2, pm (48hr post-seed)**

1) Determine in-line titer:

Read out visualizable titer (e.g. GFP), if applicable

**OR**

Pool infected cells from T175, determine cell concentration, seed  $\frac{1}{2}$  number of cells determined in **Part 1 Step 1** in 6 wells and start drug selection for drug-based titer (see **Part 2 Step 2**)

2) Seed the remaining pooled infected cells into flasks for screening and start drug selection.